



# BIOFELT

## Pledgets

Biofelt pledgets are devices made of Polytetrafluoroethylene (PTFE).

Available in two different type, **soft or hard**, they are used for the **reinforcement of sutures and to support the load**. On host tissue, which could otherwise be damaged.

Biofelt pledgets are steril, for single use, implantable and produced in different shapes and sizes with soft or hard consistency, to deal with different procedural requirements. **Custom sizes are available on request.**

**CERTIFICATION**  
CE 0373 class III

CODE	SPECIFICS	PCS. x PACK	PACKS x BOX
<b>FP0303S</b>	3x3 soft	6	36
<b>FP0306H</b>	3x6 hard	6	36
<b>FP0306S</b>	3x6 soft	6	36
<b>FP0306K</b>	3x6 hard w/holes	2	24
<b>FP0409S</b>	4,5x9 soft	6	36
<b>FP0676S</b>	6x76 soft	1	24
<b>FP2510S</b>	25x10 soft	1	24
<b>FP2525H</b>	25x25 hard	1	24
<b>FP2525S</b>	25x25 soft	1	24
<b>FP5151S</b>	51x51 soft	1	24
<b>FP1201S</b>	12x102 soft	1	12

CODE	SPECIFICS	PCS. x PACK	PACKS x BOX
<b>FP5101S</b>	5x102 soft	1	12
<b>FP0676S</b>	6x76 soft	1	12
<b>FP0615S</b>	6x152 soft	1	12
<b>FP0815S</b>	8x152 soft	1	12
<b>FP1001S</b>	102x102 soft	1	12
<b>FP1215S</b>	12x152 soft	1	12
<b>FP2515S</b>	25x152 soft	1	12
<b>FP1015S</b>	100x152 soft	1	6
<b>FP1515S</b>	152x152 soft	1	6



BIOMED® is a Bio Sud Medical Systems s.r.l. registered trademark

**BIO SUD MEDICAL SYSTEMS s.r.l.**

Via dei Fabbri, 23/25 | 70026 Modugno (BA) | ITALY

**Tel: +39 080 5355669**

**[www.biomed surgical.eu](http://www.biomed surgical.eu)**  
**[info@biomed surgical.eu](mailto:info@biomed surgical.eu)**